

Climatic Chambers.

CB-CS Serie


FDM Environment Makers Products

Climatic and Thermostatic Chambers Growth Chambers Environmental Chambers Ultrafreezers Custom Products

F.lli Della Marca S.r.I. Viale Arcangelo Ghisleri 00176 Rome (Italy) (+39) 06 29 80 42 info@dellamarca.it www.dellamarca.it

Technical specifications are subject to change.

Sizes



Advantages

Volumes from 140 up to 2300 liters

Adjustable temperature -25°C...+70°C*

Humidity range 10...98%.

Internal shelves adjustable in height

Programmable controller revoFACE or stillFACE for stable testing

Internal forced and homogeneous ventilation

High and low temperature and humidity alarms

Electrical power 220-240V /50Hz or 110V

Compliant to the ICH Q1A pharmaceutical stability test

Configuration examples



C140BXPRO PV101-140



C700BXPRO **RT100**

Description

The CB-CS series is the standard FDM climate chamber category: you can select different temperature, humidity and volume parameters.

You can also choose between a cyclic or a constant climatic chamber: the cyclic climate chamber was developed for programming different test cycles with different temperatures during the same test.

The climatic chamber is an extremely important laboratory tool also and above all because the various tests to be carried out are regulated by different standards, established by official research bodies.

The range of products to be tested is various, below are some of the most important application areas:

- cosmetic and pharmaceutical
- . food
- textile
- military/aerospace . .
- electronic components
- **Automotive**
- construction and building materials

For some categories, which require more strict and accurate parameters and also the reproduction of more "extreme" tests, we mostly refer to an environmental chamber.





The stillFACE controller is designed for stability testing.

The revoFACE programmer is designed for cycle testing.

*Except for chamber with 140 liters volume starting from -20°C up to +70°C



C700BXPRO PV101 + RT100



C1500SXPRO



Positive Range (S)		C140S	C370S	C440S	C700S	C900S	C1200S	C1500S	C2300S	
Negative Range (B)		C140B	C370B	C440B	C700B	C900B	C1200B	C1500B	C2300B	
	(S)	0+70								
Temperature range without humidity [°C]	(B)		-25+70							
Temperature performance										
Temperature uniformity depending on setpoint [± K]					0.5.	.2.5				
Temperature fluctuation depending on setpoint [± K]	0.10.5									
Average heating-up rate according to EN 60068-3-5 [K/min]					+	2				
Average cooling down rate according to EN 60068-3-5 [K/min]		-1 -0.6 -1 -0.6).6	-1				
Performance Data Climate		,								
Temperature range [°C] humidity is on					+10/	′+70				
Humidity range in according to Chart*1 [% r.h.]					10/	98				
Humidity fluctuation depending on setpoint [± %]	≤ 3									
Humidification system	Ultrasonic humidifcation system									
Water filter	Water softening filter and replaceable cartridge equipped*2									
Max Water hardness permitted		600 ppm CaCO ₃								
Water inlet	Pressure 0.2÷5 BAR Temperature 10÷40°C Pipe connection 3/4"									
Outer dimensions										
Width net [mm]		530	600	600	750	1200	1500	1500	2250	
Depth net [mm]		810	768	868	980	868	860	980	980	
Height net [mm]		950	1920	1920	2100	1920	2100	2100	2100	
Doors										
Unit doors		1	1	1	1	2	2	2	3	
Internal Dimensions										
Width [mm]		540	480	580	670	580	570	670	670	
Depth [mm]		520.5	1257	1257	1340	1257	1350	1340	1340	
Height [mm]		520.5	1257	1257	1340	1257	1350	1340	1340	
Measures										
Interior net volume [L]		130	310	370	530	810	1020	1180	1830	
Net weight of the unit (empty) [kg]		78	98	108	165	195	210	213	365	
Permitted load per rack [kg]		30	30	30	30	30	30	30	30	
Fixtures										
Number of shelves (std./max.)		3/5	3/5	3/5	3/6	6/10	6/12	6/12	9/14	
Shelf Width [mm]		400	460	460	530	460	530	530	530	
Shelf Depth [mm]		500	470	570	650	570	550	650	650	
stillFACE Controller										
Controller		Constant monitoring temperature and humidity controllers								
Interface port					Optional	: RS485				
revoFACE Programmer										
Set-up Display			Programme	e status, set u	o temperature	and humidity	, date, time ar	nd language		
Programmer			10 prog	rams and 50 s	egments each	n, adjustable f	irom 1 min. to	999 hrs.		
Calibration	Possibility to calibrate all parameters									
Interface port			Option	al: RS485 - O	ptional: Etherr	net - Wi-Fi - D	edicated mob	ile app		
Environment-specific data										
Sound-pressure level [dB(A)]					6	0				

*2The cartridge should be change every 4÷6 months at most.

All technical data are specified for units with standard equipment at an ambient temperature of +24 °C and a voltage fluctuation of ±10 %.

Positive Range (S)		C140S	C370S	C440S	C700S	C900S	C1200S	C1500S	C2300S			
Negative Range (B)		C140B	C370B	C440B	C700B	C900B	C1200B	C1500B	C2300B			
Structure and insulation												
Exterior material		Stainless steel AISI 304										
Internal material		Stainless steel AISI 304										
Insulation		CFC and HCFC free										
Shelf grill		Removable and height adjustable stainless steel plastic coated										
Ventilation												
Fan		forced ventilation										
Safety												
Temperature		Class 2 independent adjustable temperature safety devices (DIN 12880)										
Alarm		Audio-Visual										
Electrical data												
Rated Voltage [V]		220/240										
Power frequency [Hz]		50										
Nominal power [kW]	(S)	1.4	1.5	1.4	1.8	2.2	2.4	2.6	4.1			
	(B)	1.6	1.7	1.7	2.3	2.6	2.7	2.9	4.1			
Unit fuse [A]	16 2								25			
Phase (nominal voltage)	Phase (nominal voltage) 1 ~											

Main Options

- ETH100 Ethernet interface for remote connection
- RS485 Interface port RS485
- PG50200 Upgrade. 50 program 200 step for revoFACE
- PE100 Interior socket coverable universal 220-230 V
- SPS100 Interior compartment glass doors
- SPS07 Flooding system for EN1348-2007
- BDS020 Relative humidity decreasing up to 2% @ 20°C
- BSU10 Humidity extension range below the 10°C
- TE110 Power supply 110 V 50-60 Hz
- CO200 CO2 Regulation System for EN 12390-10
- ICHQ1A Equipment for pharmaceutical stability testing ICH Q1A



FD100 Cable hole installed on the side of the chamber that allows the passage of cables inside without loss of performance. Rubber cap included.





GFDM



PV101*

Double tempered

complete visibility

of the environment and stability of the

glass door for a



GP100R The reinforced stainless steel shelf differs from the standard one (pictured above) in that it has a higher supported load capacity.



WTO16 System of supply demineralizzate water tank 16 liters.




All technical data are specified for units with standard equipment at an ambient temperature of +24 °C and a voltage fluctuation of ± 10 %.

PV101' - The adoption of this accessory reduces the maximum humidity range to 80% and the minimum temperature range to -20°C.

Click and go to Online Configurator \bowtie www.dellamarca.it

Operation flowchar



B= With BSU accessory (for a limited time) D= Possible condensation inside the test chamber.

A= Standard

*1Operation Temperature Humidity Limits